

Prebid online meeting questions and reply
Tender Ref no STREAM/EQ/010/01/2025
Differential Scanning Calorimeter Machine

SI no.	Tender clause-Annexure -B	Question	Reply
1	1	Most of the application range can be covered without LN2 cooler in a temperature range from -85 degC to 500 degC	Refer corrigendum 3
2	2	Microbalance will not be required for DSC	Refer corrigendum 3
3	4	Accuracy and precision to be defined for a good specification	Refer corrigendum 3
4	5	Baseline accuracy to be defined to view weak transitions	Refer corrigendum 3
5	2	Thermocouple material is to be defined for long life apart from the type	Refer corrigendum 3
6	5	Mass loss studies and isothermal conditions will be a part of TGA/DTA equipment and not under DSC	Refer corrigendum 3
7	6	Can autosampler positions to be increased to 35 positions instead of 20	Based on current requirements, 20 position autosampler is sufficient and based on quantity.

8	6	Traceable standards of 6 will be very high and add to the cost when 2-3 is sufficient	Refer corrigendum 3
9	6	100 Alumina pans will be less and will be readily consumed	Refer corrigendum 3
10		Training duration to be looked at	Refer Corrigendum 2